

Holter Monitor



25 YEARS IN
CARDIOLOGY

www.norav.com

Holter ECG NH – 301

“Holter monitoring has never been easier”

Norav's true 12-3 lead holter system features:

- Atrial Fibrillation Detection
- 14 Day Recording
- Outstanding on-screen graphics
- Cloud Analysis
- Arrhythmia, ST, HRV & QT analysis, comprehensive report printing
- Interface with electronic health records



Additional Features:

- Review ECG, save to hard drive . Export reports
- Review and print trend
- Pacemaker detection
- Recorder with a large LCD display for setting up and ECG data preview
- Automatic & manual editing of patient data
- Network interfacing capability
- Link to ECG management software (NEMS)
- DICOM Interoperability
- Full Disclosure with easy page by page scrolling
- Wireless hook up check on PC *3
- Wireless transfer of patient info PC from to the recorder *3
- Analysis is completed and ready to be printed within minutes after downloading
- Maximum patient comfort thanks to the small Holter ECG recorder.
- Comprehensive, reliable Holter report

	NR Options			
Modules	NR-302	NR-314	NR-1207	NR-1207-3
Dimensions: 92 x 75 x 23 mm	✓	✓	✓	✓
Weight with one battery: 126 g	✓	✓	✓	✓
Battery type: R06 / AA	✓	✓	✓	✓
Operating time for long-term ECG recordings:	48H	336H *1	336H *1 / 168H *2	
ECG channels (Patient ECG Cable):				
3 channels (3 lead wires)	✓	✓	✓	✓
3 channels (5 lead wires)	✓	✓	✓	✓
3 channels (7 lead wires)	✓	✓	✓	✓
12 channels (10 lead wires)			✓	✓
ECG Analog bandwidth:	0.05 -65 Hz	0.05 -260 Hz		
ECG A2D bits: acquisition 24b/storage 12b	✓	✓	✓	✓
ECG Sample Rate:	250	250,500,1000		
Lead OFF detection:	✓	✓	✓	✓
Pacemaker detection:	✓	✓	✓	✓
Respiration (thorax impedance variations):		✓	✓	✓
Patient Physical Activity Monitoring:		✓	✓	✓
Memory capacity:	512MB	2GB	2GB	2GB
Data transfer mode:				
Flashcard reader	✓	✓	✓	✓
USB 2.0	✓	✓	✓	✓
Voice Note (Audible Diary)		✓	✓	✓
Wireless transmission of online ECG: Bluetooth® Class 1		✓	✓	✓
Simultaneous wireless transmission of online ECG+Holter recording				✓

*1 14 days using 3 channels, 250 s/sec. & a lithium battery


*2 7 days & 12 channels using a lithium battery & 250s/sec

*3 Except with NR-302

Specifications are subject to change

www.norav.com info@norav.com

Cloud Holter Scanning Center



See also our other product lines:

Resting 12 Lead ECG

Stress ECG

Telemetry ECG+SpO2*

